

Scottish Sheep and Goats Movement Document

Sheep & Goat (Identification & Traceability) (Scotland) Regulations 2006 (as amended)

This form is in triplicate. Please insert a hard surface between the form in use and the rest of the pad. Use a black ballpoint pen and write firmly to allow the carbon paper to duplicate the information.

Departure location

Section 1 of this document must be completed by the keeper moving the animals off. The departure CPH must be the CPH relating to the land the animals are moving from. You must supply your home address. **All fields in grey must be completed.** If you are moving your own animals please enter "self" in the Haulage Company field. If you do not have a Permit Number enter "None". (You must register with Animal Health to obtain a Permit Number if you are hauling stock over 65 kilometres, *approx. 40 miles*, in distance).

Please keep the pink form for 3 years and send the white & yellow copies with the animals. If you are directly exporting animals outwith Great Britain send the white copy to SAMU within 3 days of the movement.

Receiving location

Section 2 of this document must be completed by the keeper receiving the animals. **All fields in grey must be completed.** Send the white movement slip to SAMU within 3 days of the move, and keep the yellow copy for 3 years.

You may notify the movement by returning the white slip or by faxing or emailing a scanned copy of the document to SAMU. If you receive animals from a Scottish market they will notify the "on" movement for you.

Points to note

Animals arriving at an Abattoir not correctly identified may be refused entry into the food chain.

Animals losing their identity must be fitted with a red replacement eartag. No other type of tag can use this colour. Red tagged animals cannot be "live" exported.

A farm assurance number/sticker may be added to the movement document provided it does not interfere with required information.

Animals moving to and returning from a showground will require a **separate movement document** to record each movement.